

# DynaMarker, Prestain Marker For Small RNA

## Features

### ■ Prestained !

"DynaMarker Prestain Marker for Small RNA Plus" is highly visible indicator with dual colors of blue and red.

### ■ From 20 to 100 base !

The apparent sizes of the bands in the marker are 20, 30, 40, 50, 75 and 100 bases in length.


### ■ Ready to use !

It is ready-to-use mixture for acrylamide gel.


**Do not require heat treatment or denaturing agents.**

### ■ Visible during electrophoresis !

The visible prestained bands during electrophoresis enable you to run RNA electrophoresis under the optimized condition.



Electrophoresis profile of DynaMarker Prestain Marker for Small RNA Plus (5  $\mu$ L) on 10 % acrylamide containing 7.5 M urea gel with 1 $\times$  TBE buffer as running buffer.



Left: Electrophoresis profile 10 % acrylamide containing 7.5 M urea gel with 1 $\times$  TBE buffer as running buffer.  
Right: Nylon membrane after transfer

## Product Information

[Supplier : BDL]

Product Name	Size	Code	Storage
DynaMarker, Prestain Marker for Small RNA Plus	75 $\mu$ L	DM253S	-20 $^{\circ}$ C
	150 $\mu$ L	DM253	-20 $^{\circ}$ C

# DynaMarker, Prestain Marker For RNA High

## Features

### ■ Prestained !

"DynaMarker, Prestain Marker For RNA High" is highly visible indicator with dual colors of blue and purple.

### ■ For high base RNA !


The apparent sizes of the prestained bands are **200, 500, 1,000, 2,000, 4,000 and 8,000** bases in length.

### ■ Ready to use !


It is ready-to-use mixture for agarose gel. **Do not require heat treatment or denaturing agents.**

### ■ Visible during electrophoresis !

The visible prestained bands during electrophoresis enable you to run RNA electrophoresis under the optimized condition.



Electrophoresis profile 0.8% agarose - 2.2 M formaldehyde gel / 1 x MOPS buffer as running buffer.



Left: Electrophoresis profile of 0.8% agarose - 2.2 M formaldehyde gel / 1 x MOPS buffer as running buffer

Right: Nylon membrane after transfer

## Product Information

[Supplier : BDL]

Product Name	Size	Code	Storage
DynaMarker, Prestain Marker for RNA High	90 µL	DM260S	-80°C
	180 µL	DM260	-80°C

NOTE \* All products on here are for research use only, not for diagnostic use. \* All brand and product names are registered trademarks of their respective companies.  
\* Products might be changed for improvement without prior notice. \* Please contact your local distributors for orders, quote request and any inquiry.

Your Local Distributor



**PELOBiotech GmbH**

Address: Am Klopferspitz 19  
82152 Planegg / Martinsried Germany  
Phone : +49(0)89-517 286 59-0  
Fax : +49(0)89-517 286 59-88  
Email : info@pelobiotech.com



**funakoshi Co., Ltd.**

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,  
Tokyo 113-0033 JAPAN  
Phone: +81-3-5684-6296  
Fax: +81-3-5684-6297  
Email: export@funakoshi.co.jp