









## MicroStem SlideHolder™

The MicroStem SlideHolder™ enables MicroMatrix™ users to acquire images on any plate based system that is commonly used for immuno-fluorescent image acquisition.

Currently, MicroStem is working with multiple instrument providers to facilitate imaging of MicroStem products. For more information on instruments, please contact tech support at MicroStem, Inc.

### CONTACT INFO

MicroStem Inc.

10575 Roselle Street

San Diego, CA 92121

P (800) 819 2534

F (800) 819 1514

info@microstem.com

[www.microstem.com](http://www.microstem.com)

that supported MCF-7 attachment for compound toxicity assessment. The substrate identified maintains desirable adherence and morphology of MCF-7 cells for proper imaging. Cells remained immobilized during Taxol dosing at various concentrations, allowing researchers to monitor early and late stage toxic endpoints. MicroMatrix™ extracellular matrix array allows researchers to rapidly identify ECMs that influence cellular fate and function. The system uses minimal amount of samples and can be incorporated into a high throughput cell based assay format using high content imaging or other image based detection methods. Here we demonstrate the utility of MicroMatrix™ in a cellular assay when combined with a Cellomics™ vTI and its companion software. The MicroMatrix™ product offerings are a robust cell based assay tool that can be used in many applications including primary cell, stem cell and immortalized cell line assay optimization and screening.

### References:

1. Wendt MK, Allington TM, and Schiemann WP. Mechanisms of the epithelial-mesenchymal transition by TGF-beta. *Future Oncol.* 2009 Oct;5(8):1145-68. PMID:19852727
2. Collagen I Promotes Epithelial-to-Mesenchymal Transition in lung cancer cells via transforming growth factor- $\beta$  signaling. *Am J Respir Cell Mol Biol.* 2008 Jan;38(1):95-104. PMID:17673689
3. Brafman D., Shah KD, Fellner T, Chien S, Willert K. Defining Long-Term Maintenance Conditions of Human Embryonic Stem Cells With Arrayed Cellular Microenvironment Technology. *Stem Cells Dev.*, 18(8): 1141-1154, 2009 PMID: 19327010
4. Flaim, CJ., Teng, D., Chien, S. and Bhatia, SN. Combinatorial signaling microenvironments for studying stem cell fate. *Stem Cells Dev* 17:29-39, 2008, PMID: 15782209



Product Name	Catalog
MicroMatrix™ 36	MM-012011
MicroMatrix™ 96	MM-022011
MicroMatrix™ 192	MM-032011
SlideHolder™	SH-042011