

Cancer Related Inhibitors

For more information : <http://www.funakoshi.co.jp/exports>

Inhibitor for DNA Helicase


Heliquinomycin

Heliquinomycin inhibits DNA helicase from HeLa cell in a non-competitive manner with the inhibition constant (Ki) of 6.8 mM. The topoisomerase II and I enzymes are inhibited at 30 µg/ml and 100 µg/ml of heliquinomycin, respectively. Heliquinomycin inhibits the growth of HeLa S3, KB, LS180, K562 and HL60 human tumor cell lines at IC₅₀ values of 0.96 to 2.8 µg/ml. Heliquinomycin inhibits both DNA and RNA synthesis in cell culture but does not inhibit protein synthesis. HeLa S3 cells are arrested at the G2/M phase by heliquinomycin.

This product is licensed under JP patent NO.3490154

Reference

- Chino, M., *et al.*, *J. Antibiotics*, 49, 752~757 (1996).
Chino, M., *et al.*, *J. Antibiotics*, 50, 143~146 (1997).
Chino, M., *et al.*, *J. Antibiotics*, 50, 781~784 (1997).
Chino, M., *et al.*, *J. Antibiotics*, 51, 480~486 (1998).




[Manufacturer : IMC]

Product Name	Size	Catalog #	Storage
Heliquinomycin	0.2 mg	10665	4°C

Inhibitor for Heparanase

Heparastatin <SF4>

Heparastatin (SF4) inhibits recombinant human heparanase from human melanoma A375M cells transfected with pBK-CMV expression vectors containing the heparanase cDNA with IC₅₀ 1.02 µM. Heparastatin (SF4) inhibits β-D-glucuronidase from bovine liver with IC₅₀ 6.5 x 10⁻² µM. Heparastatin (SF4) (100 µM) completely inhibits the enzyme activity of recombinant heparanase of murine mammary epithelial cells (NMuMG) transfected with a mouse heparanase expression vector pcDNA3.1(-)-Hygro-Hep at 0.15 µg/ml in a *in vitro* HS degradation assay. Heparastatin (SF4) inhibits heparan sulfate (HS) chain degradation of HSPGs of Matrigel by heparanase of the LPS-treated microglial lysates from the forebrain cells of Wistar rats and the *in vitro* transmigration of microglia through the Matrigel-coated insert in a dose-dependent manner. Heparastatin (SF4) markedly inhibits degradation of HS by heparanase in the nucleus translocated from the cytoplasm of the calcium-induced human esophageal keratinocyte cells and keratinocyte differentiation at 100µM. Heparastatin (SF4) markedly inhibits in a dose-dependent manner experimentally induced pulmonary metastasis of the B16BL6 in mice. Inhibitory ratio by *ex vivo* treatment with 50 µg/ml of Heparastatin (SF4) is 90.8%. Heparastatin (SF4) shows 57% inhibition of lung metastasis of 3LL cells by s.c. inoculation in mice with i.v. administration of 100 mg/kg/day for 5 days.



[Manufacturer : IMC]

Product Name	Size	Catalog #	Storage
Heparastatin <SF4>	0.2 mg	11829	4°C

Inhibitor for Protein Phosphatase 2A


Cytostatin

Cytostatin inhibits the adhesion of B16 melanoma cells to laminin and collagen type IV in a dose dependent manner but not to fibronectin. The IC_{50} values are 1.3 $\mu\text{g/ml}$ to laminin and 1.4 $\mu\text{g/ml}$ to collagen. The administration of cytostatin inhibits metastases of B16-F10 markedly. The inhibitory ration is about 60 to 70% at 1.25 mg/kg/day. Cytostatin inhibits protein phosphatase 2A with an IC_{50} of 0.09 $\mu\text{g/ml}$ in a non competitive manner against a substrate, p-nitrophenyl phosphate, but it has no apparent effect on other protein phosphatases including protein phosphatase 1, protein phosphatase 2B and alkaline phosphatase even at 100 $\mu\text{g/ml}$.

This product is licensed under JP patent NO.367455

Reference

- Amemiya, M., *et al.*, *J. Antibiotics*, 47, 536~540 (1994).
 Masuda, T., *et al.*, *J. Antibiotics*, 48, 528~529 (1995).
 Kawada, M., *et al.*, *Biochim. Biophys. Acta.*, 1452, 209~217 (1999).
 Kawada, M., *et al.*, *Int. Immunopharmacol.*, 3, 179~188 (2003).



[Manufacturer : IMC]


Product Name	Size	Catalog #	Storage
Cytostatin	0.2 mg	10664	-20°C

Rubratoxin A

Rubratoxin A is a mycotoxin produced by *Penicillium*. While this compound inhibits protein phosphatase 2A (PP2A) in a competitive manner with the inhibition constant (K_i) of 28.7 nM (IC_{50} = 6.1-15.1 nM), it does not inhibit protein phosphatase 1, protein phosphatase 2B, protein tyrosine phosphatase 1B, alkaline phosphatase, and calf intestine phosphatase(CIP) up to 200 μM . Rubratoxin A shows cytotoxicity in cultured mammalian cells with IC_{50} values of around 15 μM and induces overphosphorylation of PP2A substrate proteins.

Reference

- Wada, S., *et al.*, *Cancer Sci.*, 101, 743~750 (2010).



[Manufacturer : IMC]

Product Name	Size	Catalog #	Storage
Rubratoxin A	0.2 mg	10663	4°C

NOTE

※ All products here are research use only, not for diagnostic use.
 ※ Specs might be changed for improvement without notice.

※ Company name and product name are trademark or registered mark.
 ※ Please contact your local distributors for orders, quote request and inquiry.

Your Local Distributor



PELOBiotech GmbH

Address: Am Klopferspitz 19
 82152 Planegg / Martinsried Germany
 Phone : +49(0)89-517 286 59-0
 Fax : +49(0)89-517 286 59-88
 Email : info@pelobiotech.com



funakoshi Co., Ltd.

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,
 Tokyo 113-0033 JAPAN
 Phone : +81-3-5684-6296
 Fax : +81-3-5684-6297
 Email : export@funakoshi.co.jp