

Cell Encapsulation Protocol

Under aseptic conditions add 9 mL of sterile deionized water to *mimsys*[®] G vial, dissolve it under agitation in a 37°C water bath. **Work with 3D cell culture material at physiological temperature!**

Prepare your cell suspension in 1 ml of your favorite culture media.

Great advantage when working with a 3D cell culture system without adding cytotoxic crosslinks!

Mix the 1 ml cell suspension with the *mimsys*[®] G, and transfer it to a mold or directly to your well plate. **You can remove the mold in 5 minutes!**

Cover the encapsulated cell on *mimsys*[®] G hydrogel with culture media and incubated your experiment! **Change the culture media with your normally routine during all experience duration.**

**Culture media for crosslink!
Non-cytotoxicity, full biocompatibility!**

More Protocols Available

- Viability assay by MTS
- Live/Dead assay and cell morphology imaging by fluorescent microscopy
- Standard histology staining protocols: H&E and chondrogenic differentiation

Intended Use

Research Use Only

Shipment


At room temperature

Storage

At 15°C to 20°C
Keep in a dry place or store in inert atmosphere is recommended

version 1.1 updated on January 2014 © All Rights Reserved

irisbiosciences is a trademark of Stematters, Biotecnologia e Medicina Regenerativa S.A.



3D Cell Culture

mimsys[®] G

xeno-free & nutrient permeable hydrogel for 3D cell culture

mimsys[®] G is a xeno-free, non-immunogenic, easy to handle hydrogel used for cell encapsulation in 3D experiments *in vitro* or *in vivo*. *mimsys*[®] G is nutrient permeable and allows long term 3D cell culture, with a viability of at least 21 days of proliferation

The transparency of *mimsys*[®] G hydrogel makes it suitable for established laboratorial cell assays, such as cell differentiation, microscopy evaluations, viability protocols and histological procedures

Characteristics

Xeno-free

Transparent

Water soluble

Defined composition: methacrylated gellan gum is based on gellan gum which is extensively used in food industry providing no-cytotoxicity and no-immunogenic response

Gelation in less than 5 min by ionic crosslink with cationic solution (such as media or PBS)

Gelation by photo-initiation with UV light for higher spacial and temporal control

Easy to handle at room temperature and at 37°C

Allows oxygen and nutrient flow through the hydrogel up to the cells

Benefits

mimsys®G accurately mimics the 3D cell environment:

- Cell-to-cell communication is improved
- Cell morphology is closer to its native one
- Right flow of nutrients through the hydrogel

- No necrosis: oxygen can naturally penetrate the gel
- High reproducibility of prepared hydrogels due to the robustness and easy handling of the product
- Compatible with stem cells and cell lines 3D culture for long term culture
- Easy to use in *in vivo* experiments

Appropriate for established downstream applications:


- Optical microscopy observation of your 3D cell culture
- Viability assay, such as MTS
- Fluorescence microscopy for Live/Dead and morphology assays
- Histology assays of hydrogels from *in vivo* studies

Other application range:


Drug delivery, diagnostic marker delivery and surface coatings

Save time in your research:
Fast gelation time – 5 min

Save money in your experiment:
1 vial suitable for 2x96 well plates



Cell morphology in long term culture (21 days)



High cell viability



Customize 3D hydrogels with mimsys®G

Customize your experiment!

Hydrogel with the shape and volume that you want (well-plates, petri-dish, etc)

Easy supplementation with growth factors, matrix proteins, or other molecules specific for your work

Injectable hydrogel with *in situ* gelation to improve the control in your *in vivo* studies



Related products

060230 · hSVF
human Stromal Vascular Fraction

060231 · hASC
human Adipose derived Stem/
Stromal Cells

060232 · hASC xeno-free
human Adipose derived Stem/
Stromal Cells – xeno-free

060233 · hSVF xeno-free
human Stromal Vascular Fraction –
xeno-free