


PP:Autoclavable
 PE:Non-autoclavable
 POM:Non-autoclavable
 Nitrile rubber:Autoclavable
 Zirconia beads:Autoclavable




DISPOSABLE HOMOGENIZER BioMasher®


Made in Japan

Distributed by PELOBIOTECH in Austria, Germany and Switzerland - call us +49 89 517 286 59 0 or send an Email: info@pelobiotech.com

BioMasher-I




3 The pestle is necessary.




10,000~15,000xg 10~30sec


BioMasher-II



BioMasher-III




3 The pestle is not necessary.



6,000~10,000xg 10~30sec

BioMasher-IV



	BioMasher-I	BioMasher-II	BioMasher-III	BioMasher-IV
Description	For many samples	For increase of the collection amount	For increase of the collection amount For many samples	For hard samples For many samples
Collection amount	Small	Large	Large	Large
Heat	Heatless	Controllable	Controllable	Heat
Device	Centrifuge	Unnecessary	Centrifuge	Cell crusher
Soft samples	⊙	⊙	⊙	○
Hard samples	△	⊙	⊙	⊙
Yeast	×	×	×	⊙
Cultured cell	△※	×	△	○

※Please inquire for the usage.