

# Human Male AB Serum (Off-the-clot)



Distributed by



**PELO**Biotech

More than Competence<sup>4</sup> Cells

P +49 (0)89 517 286 59 0  
www.pelobiotech.com

## HUMAN AB SERUM: TRUE OFF-THE-CLOT QUALITY TESTED SERUM FOR YOUR NEEDS

Akron's AB serum pools are produced from true off-the-clot (OTC) human serum, from healthy males who test negative for HBsAg, HCV and HIV 1/2 and is ideal for culturing cells associated with the human immune system.

### Features

- Male AB serum is low in estrogen and HCG
- Viral tested and found negative for HBsAg, HCV and HIV 1/2.
- Sterility testing is performed according to USP guidelines.
- Gamma irradiation and heat inactivation available upon request



### Product

### Cat#

Human AB Serum (True, off-the-clot)	AK9905
Human AB-Serum (Converted, off-the-clot)	AK9340

### Related products

- Platelet Lysate
- Human Serum
- Cord Blood Serum
- Affinity-stripped Plasma



sales@akronbiotech.com  
www.akronbiotech.com

**AKRON**<sup>TM</sup>  
BIOTECH